




**Instructions for adding the Osteoscan Template to Best Practice  
(tested on Saffron Edition Build No: 1.11.0.924).**

1. Please download the template file (**OsteoscanBP.rtf**) from our website, place the downloaded file on the **“Desktop”**


Alternatively download the Zip file, open the Zip file and extract the template from inside, place the extracted template on the Desktop.

2. Open Best Practice, login as a user with high privileges (ie: Practice Manager)


3. Open the **“Word Processor”** (it’s in the upper menu, Utilities then Word Processor)




4. Once inside the **“Word Processor”**, go the upper menu, click on Templates and then **“Import Template”**



5. On the window that pops up titled **“Open”**, select the **“Desktop”** from left sidebar, then select the file we downloaded in step 1, and click **“Open”**



6. Within the Word Processor, on the top menu, select File then **“Save as”**, a new window will pop up asking for the Template Name, type in **“Osteoscan”** and tick the box that says **“Available to all users”**



7. Hit **“Save”** and close all Windows, now the template is available for anyone to use